

T34™

Deliver Syringe Infusion Anywhere

Provide complete small-dose infusion therapy safely and economically wherever your patients receive care.

The T34™ ambulatory syringe pump from CME contains many features of larger syringe pumps in one highly compact, lightweight design.

The T34™ pump is easy to use for inpatient and outpatient applications.



Applications

- Medications/analgesics
- Palliative care
- Brittle Asthma
- Chemotherapy
- General fluid replacement

Intended Use

- Intravenous
- Subcutaneous
- Percutaneous
- Epidural (in close proximity to nerves)
- Intra-operative site
- Intra-arterial



Features & Benefits

Automatic Syringe Detection

Automatically detects any syringe brand using 3 independent sensor system.

Unique algorithm detects your drugs volume to reduce dosing errors.

Fully Programmable Infusion

Deliver your preferred continuous or intermittent protocol.

Programmable Flow Rates

0.1 – 650 ml/h range to meet your every needs.

Programmable Infusion Time

Instead of rate, you can select from 10 minutes to 100 hours of default infusion time. The T34™ syringe pump then automatically calculates flow rate.

Configure Group Settings

BodyComm™ software with serial port allows you to customize an infusion protocol and easily program multiple T34™ pumps with the same settings.

Inpatient/Outpatient Accessories

A lockbox is available for pole-mounted use, and/or a carrying pouch for home care.

BodyComm™ Software

Review patient history with charts and 1000 events log or manage pump settings on your PC.

Air & Mobile Medicine

CME infusion pumps provide a complete range of infusion onboard any air or ground vehicle. CME pumps are RTCA-DO160 compliant.

True Portability

Size

167 x 66 x 39 mm.

The T34™ syringe pump fits anywhere.

Weight

277 g

Patients can carry it wherever they go.

Battery-Powered

Operates using one 9V Alkaline battery.

The battery lasts for 7 full devilsries.

Flexible Configuration

Supports Multiple Languages

Over 20 languages available to meet your needs and preferences.

Safety & Security

Full Alarm Protection

Helps prevent events and instill confidence in patients.

Descriptive Status/Alarm Display

Easy-to-read backlit LCD.

Easy-Scroll Event Review

Press INFO key for hour-by-hour event log.

3-Level Access Codes

Prevent accidental or unauthorized program changes.

